
- Oscillateur à diode GUNN
- GUNN diode oscillator
- GUNN-Oszillator
- Oscillatore a diodo GUNN
- Oscilador de diodo GUNN

ORITEL OSG 100

FRANÇAIS
ENGLISH
DEUTSCH
ITALIANO
ESPAÑOL

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English

Meaning of the symbol 

CAUTION! Consult the operating manual before using the instrument.

In this operating manual, the instructions preceded by this symbol may cause personal injury or damage the instrument and installations if they are not scrupulously abided by or carried out.

This instrument is sensitive to electrostatic discharges.

You have just bought an **ORITEL OSG 100 Gunn diode oscillator** and we thank you for the confidence you have in our products.

In order to get the best use out of your instrument:

- carefully **read** this operating manual,
- **follow** the precautions for use.

PRECAUTIONS FOR USE

- Because of the diode's sensitivity to electrostatic discharges, considerable precautions should be taken when handling this oscillator.
- When the Gunn oscillator is connected to the power supply, any inversion of the power supply output connectors may lead to the Gunn diode being damaged (see § 3.3).
- Comply with the Gunn oscillator's supply voltage: 10 V_{dc}
- The oscillator must be fitted along with an ORITEL ISO 100 isolator, in order to be protected against load impedance variations which could bring about cases of frequency and level instability.
- The ferrite isolator shall be fitted in the correct direction, with the arrow pointing towards the load.
- The connecting leads and accessories shall comply with the applicable standards and shall be designed for an overvoltage category rated voltage which is at least equal to that of the circuits on which the measurements are carried out.
- This component's mass shall be earthed, as shall all the instruments and components which are connected to the bench and the AC network.
- No voltage of more than 30 V in relation to the earth shall be present on the bench.

GUARANTEE

Our guarantee is applicable for **twelve months** after the date on which the equipment is made available (extract from our General Conditions of Sale, available on request).

TO ORDER

ORITEL OSG 100 Gunn diode oscillator **P01275307**

Delivered with this operating manual

CONTENTS

1. PRESENTATION	9
2. DESCRIPTION	9
3. USE	10
3.1 Summary of the components making up the ORITEL BDH R100 basic educational bench	10
3.2 How to use the EASYFIX quick fastening lug	10
3.3 Connecting the Gunn oscillator to the ORITEL CF 204A power supply unit	10
3.4 A summary of the Gunn-effect theory	11
3.5 Determining the current-voltage characteristic by experiment	11
3.6 Determining the displacement-frequency characteristic of the Gunn oscillator by experiment	12
3.7 Determining the power-frequency characteristic of the Gunn-effect oscillator	12
4. CHARACTERISTICS	13
5. MAINTENANCE	13
5.1 Cleaning	13
5.2 Metrological inspection	13
6. APPENDIX	33
6.1 Gunn diode oscillator	33
6.2 Front panel of the ORITEL CF 204A power supply unit	34
6.3 EASYFIX quick fastening lug	34

1. PRESENTATION

The **ORITEL OSG 100 Gunn diode oscillator** is a microwave energy source specially designed for the ORITEL BDH R100 educational bench, manufactured using a rectangular waveguide in accordance with the R100/WR90 standard (covering the frequency range between 8.5 and 9.6 GHz).

It comprises a filter and a length of waveguide, closed at the end by a short circuit and on which a Gunn diode has been adapted. The Gunn diode, energised by the ORITEL CF 204A power supply unit, generates a signal whose frequency depends on the volume of the cavity formed by the waveguide. The short circuit, actuated by a precision micrometer, enables the frequency to be set and the Gunn diode oscillator's displacement/frequency response curve to be determined experimentally.

2. DESCRIPTION

See § 6.1

- ① Micrometer enabling the oscillator's tuning frequency to be set.
- ② Female BNC-type power supply connector.
- ③ Knurled tightening knob for the EASYFIX quick fastening lug.

3. USE

The manipulations described in § 3.5, 3.6 and 3.7 require at a minimum that the components mentioned on the flow chart appearing in each of the paragraphs be assembled. Generally speaking, these components are a part of the standard supplies delivered with the ORITEL BDH R100 basic educational bench.

To use these components, refer to their respective operating manuals.

3.1 Summary of the components making up the ORITEL BDH R100 basic educational bench

ORITEL OSG 100 Gunn diode oscillator	P01275307
ORITEL ISO 100 ferrite isolator	P01275308
ORITEL MOD 100 pin diode modulator	P01275309
ORITEL ATM 100 variable attenuator with micrometer adjustment	P01275310
ORITEL OND 100 curve wavemeter	P01275311
ORITEL LAF 100 measuring line	P01275312
ORITEL ADZ 100/3 three-screw impedance adaptor	P01275313
ORITEL TGN 100 waveguide-to-coax transition element	P01275314
ORITEL DEN 100 coaxial detector	P01275315
ORITEL CHG 100 adapted load	P01275316
ORITEL CC 100 short circuit plate	P01275317
ORITEL SUP 100 waveguide support	P01275318

■ Other components (R100/WR90 standard waveguide)

ORITEL AFR 100 EASYFIX quick fastening lug adaptor	P01275301
ORITEL RS 100 bench support rail	P01275303
ORITEL ST 613 thermocouple probe (for ORITEL MH 600)	P01285101
ORITEL CF 204A power supply unit for Gunn diode oscillator	Please consult us
ORITEL IR 205 SWR indicator	P01270501
ORITEL RD 100 displacement sensor for measuring line	P01275302
ORITEL ANC 100/15 (15 dB) horn antenna	P01275304
ORITEL CGX 100/20 waveguide cross coupler	P01275305
ORITEL IRIS 20-30 20 and 30 dB) coupling irises	P01275306

To use these components, refer to the corresponding operating manual.

3.2 How to use the EASYFIX quick fastening lug

See diagram § 6.3

- Arrange the two components to be assembled, taking care to correctly place the centring studs ① opposite the holes ②
- Simultaneously tighten the two knurled knobs ③

3.3 Connecting the Gunn oscillator to the ORITEL CF 204A power supply unit

Prior to energisation (*see photo § 6.2*):

1. Turn the "Gunn oscillator" adjusting knob ③ anticlockwise as far as the limit stop (zero setting).
2. Connect filter FIL 100 to the female BNC connector of the GUNN oscillator (see §6.1).
3. Connect the free end of the FIL 100's BNC power supply lead to the corresponding power supply female BNC connector ④.
4. Press the network push button ① and the indicator light ② should light up.
5. Energise the power supply unit.

6. Set the Gunn oscillator's power supply voltage to the voltage desired by means of the control knob ③. The value of this voltage is displayed on the galvanometer ⑩. The current supplied by the power supply unit is directly displayed on the galvanometer ⑨.


CAUTION!

When the Gunn oscillator is connected to the power supply, any inversion of the power supply output connectors ④ and ⑥ may lead to the Gunn diode being damaged.

3.4 A summary of the Gunn-effect theory

Gunn-effect devices behave like an electron oscillator inside a block of gallium arsenide (GaAs), formulated as below.

It can be considered as a diode, but it does not include any p-n junction and, therefore, has no rectification characteristic.

Its operation is based on the presence of a space charge which passes through the semi-conductor from the cathode to the anode during each oscillation current cycle. The effect is thus associated with a transit time.

Fitted into a resonant cavity, the Gunn-effect diode is used as an active component for producing microwaves. It can therefore be considered as a device converting a low direct voltage into microwave energy.

In that case, the oscillation frequency is determined by the volume of the cavity.

3.5 Determining the current-voltage characteristic by experiment

The Gunn diode's «Current-Voltage» response curve has a part with a negative dynamic range (see Fig. 3-1)


Figure 3-1

The assembly set up in figure 3-2 makes it possible to determine this negative dynamic range part experimentally and show up its characteristic appearance.


Figure 3-2

Once the assembly has been set up, proceed as follows:

- Turn the «OSC GUNN» voltage control knob on the CF 204 power supply unit anticlockwise as far as the limit stop (supply voltage ≥ 0 V).
- Progressively increase the voltage by increments of 0.5 V, without ever exceeding 11 V, and determine the corresponding current experimentally (direct reading on the CF 204 power supply unit's galvanometers).
- Plot the «Current-Voltage» characteristic of the Gunn effect diode and compare it to that in figure 3-1.

3.6 Determining the displacement-frequency characteristic of the Gunn oscillator by experiment

The frequency of the wave generated by the Gunn oscillator is measured so that the «displacement-frequency» characteristic can be determined; indeed, the resonance frequency depends on the volume of the cavity determined by the position of the mobile short circuit.

The frequency is set as follows:

- Consult the calibration curve of the OND 100 wavemeter and position its micrometer at the value corresponding to the frequency desired.
- Turn the Gunn oscillator's control micrometer until you obtain the frequency tuning and, consequently, the «displacement-frequency» determination.
- Carry the «displacement-frequency» correspondance points over onto a graph and plot the curve.

NOTE:

- The micrometer shall rotate clockwise if the oscillation frequency is too low, and in the opposite direction if it is too high.
- In the vicinity of the resonance, it is recommended that the micrometer be turned slowly.

3.7 Determining the power-frequency characteristic of the Gunn-effect oscillator

In the same way, the power is determined in relation to the oscillation frequency in order for its graph to be plotted.

The assembly shown in figure 3.3 is set up.


Figure 3.3

- Position the ATM 100 attenuator on 10 dB.
- Remove the square signal modulation on the CF 204 and suppress the attenuation by means of the level control knob.
- Energise the CF 204.

- Read the power indicated by the MH 600 milliwattmeter.
- Multiply this value by 10 in order to obtain the power actually supplied by the Gunn oscillator, so as to take into account the 10 dB attenuation introduced by the ATM 100 attenuator.
- Measure the frequency supplied by the Gunn oscillator using the OND 100.
- Carry the power measured in relation to the frequency over onto a graph.
- Determine several points in the 8.5 GHz - 9.6 GHz frequency band.

4. CHARACTERISTICS

Frequency: 8.5 to 9.6 GHz

Output power: > +17 dBm

Supply voltage: 10 V_{DC}

Waveguide: R100

Flange: UBR 100

Power supply connector: female BNC

Weight: 535 g

Dimensions: 53 x 59 x 225 mm

Electromagnetic compatibility:

Emissions and immunity in an industrial setting compliant with EN61326-1

5. MAINTENANCE

-  **For maintenance, use only specified spare parts. The manufacturer will not be held responsible for any accident occurring following a repair done other than by its After Sales Service or approved repairers.**

5.1 Cleaning

Before assembling the educational bench's microwave components, make sure that no dust is present inside the waveguide.

-  **CAUTION!**
Do not use a pressurised air jet for fear of damaging certain components.
Do not use solvent, and take precautions when dusting.

5.2 Metrological inspection

-  **Like all measuring or testing devices, the instrument must be checked regularly.**
 This instrument should be checked at least once a year. For checking and calibration, contact one of our accredited metrology laboratories (information and contact details available on request), at our Chauvin Arnoux subsidiary or the branch in your country.

■ Repair

For all repairs before or after expiry of warranty, please return the device to your distributor.

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ITALIA - Amra SpA
SCANDINAVIA - CA Mätssystem AB
SCHWEIZ - Chauvin Arnoux AG
UNITED KINGDOM - Chauvin Arnoux Ltd
CHINA - Shanghai Pu-Jiang - Enerdis Instruments Co. Ltd
USA - Chauvin Arnoux Inc - d.b.a AEMC Instruments