C 320

CONDUCTIMÈTRE / CONDUCTIMETER

Measure up

Français

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English

You have just purchased a **C 320 Conductimeter**, and we thank you.
For best results with your device:

- **Read** these operating instructions carefully,
- **Observe** the precautions of use.

MEANING OF THE SYMBOLS USED

	WARNING, risk of DANGER! The operator must refer to this user's manual whenever this danger symbol appears.
	The CE marking indicates conformity with European directives, in particular LVD and EMC.
	The rubbish bin with a line through it indicates that, in the European Union, the product must undergo selective disposal in compliance with Directive WEEE 2002/96/EC. This equipment must not be treated as household waste.

CONTENTS

1. TECHNICAL CHARACTERISTICS.....	28
2. DESCRIPTION OF THE INSTRUMENT.....	29
2.1 FRONT PANEL.....	29
2.2 REAR PANEL.....	31
2.3 BOTTOM.....	31
3. PRINCIPLE OF THE CONDUCTIVITY MEASUREMENT.....	31
4. USE.....	33
4.1 CONNECTION AND STARTING UP.....	33
4.2 TEMPERATURE ADJUSTMENT AND MEASUREMENT.....	33
4.2.1 <i>AUTOMATIC CORRECTION</i>	34
4.2.2 <i>MANUAL CORRECTION</i>	35
4.3 CONDUCTIVITY ADJUSTMENT AND MEASUREMENT.....	36
4.3.1 <i>CHOICE OF CONDUCTIVITY RANGE</i>	36
4.3.2 <i>ADJUSTING THE CELL COEFFICIENT</i>	37
4.3.3 <i>VIEWING THE CELL COEFFICIENT</i>	38
4.3.4 <i>ADJUSTING THE TEMPERATURE COEFFICIENT</i>	38
5. TABLE: CONDUCTIVITY OF POTASSIUM CHLORIDE SOLUTIONS.....	39
5.1 CONDUCTIVITY VS. TEMPERATURE.....	39
5.2 CONDUCTIVITY OF REFERENCE SOLUTIONS OF POTASSIUM CHLORIDE AT 20°C.....	40
5.3 ADJUSTING RECORDER OUTPUT.....	41
5.3.1 <i>FACTORY OUTPUT SETTING</i>	41
5.3.2 <i>ADJUSTING THE PARAMETERS OF THE ANALOG OUTPUT</i>	41
5.4 RS 232/V24 SERIAL DIGITAL OUTPUT.....	43
5.5 RESET.....	44
5.5.1 <i>RESET PROCEDURE:</i>	44
5.5.2 <i>VALUES AFTER RESET:</i>	44
6. IF SOMETHING GOES WRONG.....	45
7. SAFETY RULES.....	46
8. MAINTENANCE.....	47
8.1 CLEANING.....	47
8.2 METROLOGICAL CHECK.....	47
8.3 REPAIR.....	47
9. TO ORDER.....	47
9.1 DELIVERY CONDITION.....	47

1. TECHNICAL CHARACTERISTICS

MEASUREMENT RANGES		
Conductivity (6 ranges)	Resolution	RMS voltage across the terminals of the cell
0 – 200mS/cm	0,1mS/cm	80mV – 650Hz
0 – 20mS/cm	0,01mS/cm	80mV – 650Hz
0 – 2000µS/cm	1µS/cm	80mV – 650Hz
0 – 200µS/cm	0,1µS/cm	80mV – 650Hz
0 – 20µS/cm	0,01µS/cm	80mV – 90Hz
0 – 2000nS/cm	1nS/cm	80mV – 90Hz
Temperature		
Resolution		
Automatic correction (by platinum resistance probe, 100Ω at 0°C-IEC751)		
-10 – 200°C	0.1°C	
Manual correction		
0 – 200°C	0.1°C	
DISPLAY		
LCD Display unit, 2000 points, 3 ½ digits, height 17.8mm, for the conductivity measurement.		
Alphanumeric, 8 characters, height 7mm, for the temperature and messages.		
KEYPAD		
Keypad with 9 touch keys		
OUTPUT		
Recorder output 0-5V, adjustable.		
24V RS232 output for communication with PC.		
CONDITIONS OF USE		
Temperature	0 to 50°C	
Relative humidity	5 to 80% at 35°C	
Category of use	2	
Maximum altitude	2000m	

COMPLIANCE WITH STANDARDS	
Safety	With the power supply provided: EN61010-1 Normal environmental conditions (indoor use) Degree of pollution 2 (non-conducting pollution)
CEM	EN61326, class B instrument, intermittent operation.
DIMENSIONS / WEIGHT	
275mm x 208mm x 51mm / 780g	
POWER SUPPLY	
The C 320 has a 230V, 50/60Hz power supply unit	
Output of the power supply unit: 9V $\overline{\overline{=}}$	
Consumption: 0,7W	

2. DESCRIPTION OF THE INSTRUMENT

2.1 FRONT PANEL

- LCD display 2000 points
- Alphanumeric display, 8 characters
- 9 keys:

On/Off key

Cell coefficient adjustment key.

English

Temperature coefficient adjustment key.

Temperature adjustment key.

Range selection key (automatic or manual).

Analog output adjustment key.

Decrementation key.


Incrementation key.


Validation key.

2.2 REAR PANEL

2.3 BOTTOM

Identification label.

3. PRINCIPLE OF THE CONDUCTIVITY MEASUREMENT

The electrical conductance, G , of a liquid depends on the geometry of the volume of liquid involved. It is the reciprocal of the resistance, R . It is expressed in Siemens (S).

$$G = 1 / R ; 1 S = 1 \Omega^{-1}$$

The conductivity γ is the reciprocal of the resistivity ρ . It is expressed in Siemens per centimetre. The conductivity γ is a specific characteristic of the liquid. It is the conductivity of a column of liquid 1cm long having a cross-sectional area of 1cm^2 . The conductivity of an electrolyte depends only on its nature and its concentration.

$$\gamma = 1 / \rho \quad ; \quad 1 \text{ S.cm}^{-1} = 1 (\Omega.\text{cm})^{-1}$$

Just as the resistance and the resistivity are related by the formula: $R = \rho (l / S)$

The conductance and the conductivity are related by the formula: $G = \gamma (S / l)$

Where

G conductance in Siemens

γ conductivity in S.cm⁻¹.

S area in cm².

l length in cm.

The conductivity of an electrolyte can be measured simply by maintaining a fixed alternating voltage on the terminals of the cell and measuring the current that flows through it.

Characteristics of the measurement cell:

The conductometry cell comprises two square plates (5mm on a side) of platinum-plated platinum, plane, parallel, and 4mm apart. The smooth platinum plates are platinum-plated, in other words covered with platinum dendrites to form a large porous surface area. The purpose of this operation is to limit polarization of the electrodes.

The measured conductance of an electrolytic solution depends not only on the conductivity of the volume of solution exactly between the plates of the sensor, but also on the solution near this cubic volume, because the measurement cell is immersed in a larger volume of solution.

Calibrating the measurement system:

In practice, the « cell constant » of the measurement line can be determined only by calibration, by immersing the cell used in a standard solution having a known conductivity at a given temperature. (see page 39).

4. USE

4.1 CONNECTION AND STARTING UP

Unpack the instrument, the power supply unit, the calibration solution, and the operating manual.

Connect the conductivity cell (not supplied with the instrument).

Without temperature, BNC connector.

With temperature (or temperature only), lockable 5-pin plug.

Connect the power supply unit to the "power supply" connector.

(The C320 must be used with the power supply unit provided).

Press the Off/On key:

The display unit indicates conductivity measurement.

The alphanumeric display alternates between the conductivity range and the temperature.

Press the Off/On key: to switch the instrument off.

4.2 TEMPERATURE ADJUSTMENT AND MEASUREMENT

The temperature correction refers the conductivity values to 25°C.

On the C320, the temperature correction can be applied manually or automatically.

This value is indicated by the alphanumeric display.

4.2.1 AUTOMATIC CORRECTION

To use this function:

Place the probe in the liquid to be tested.

Automatic correction replaces the action of the manual corrector, which is then inoperative.

It must not be forgotten that the response of the probe is not instantaneous, and that the response time depends on the temperature gradient.

After stabilization, the temperature measurement read on the alphanumeric display may be slightly different from that indicated by a reference thermometer.

If this happens, the Pt 100 probe can be calibrated as follows:

Example: displays TEMP 22,5°C instead of 20°C

10.00 **MICRO/CM**
TEMP 22.5

Press the

key.

This instrument changes to temperature calibration mode

22.5 **CAL TEMP**

The numerical display then indicates the temperature.

Press the:

or

key

To set the desired value

20.00 **CAL TEMP**

Then on the:

key to validate.

The instrument then returns to the conductivity measurement mode.

10.00 **MICRO/CM**

TEMP 20.0

Note: maximum offset: $\pm 3^{\circ}\text{C}$

Beyond this offset, the alphanumeric display indicates:

REFUSE

4.2.2 MANUAL CORRECTION

This function can be used only if the temperature probe is not connected.
For the adjustment, the temperature of the product to be measured must be displayed.

The factory setting of this parameter is **MAN 25.0**.

Example: solution or standard at 20.0°C:

10.00 **MICRO/CM**
MAN 25.0

Press the:

key

The instrument switches to manual temperature adjustment

25.00 **REG MAN**

The numerical display then indicates the manually set temperature.

or

key

Press the:

To set to the desired value

20.00 **REG MAN**

key to validate.

Then on the:

The instrument then returns to the conductivity measurement mode.

10.00 **MICRO/CM**

MAN 20.0

Note: adjustment from: 00.0 to 200°C

4.3 CONDUCTIVITY ADJUSTMENT AND MEASUREMENT

4.3.1 CHOICE OF CONDUCTIVITY RANGE

The C 320 has 6 conductivity measurement ranges, with automatic switching or locking in one of the ranges.

To access the locked ranges or automatic switching:

Press the:

key.

The numerical display indicates:

0.00

The alphanumeric display indicates:

Range

Then, successively:

Locked range 0-200mS/cm

200 MILLI

Locked range 0-20mS/cm

20 MILLI

Locked range 0-2000 μ S/cm

2000 MICRO

Locked range 0-200 μ S/cm

200 MICRO

Locked range 0-20 μ S/cm


20 MICRO

Locked range 0-2000nS/cm

2000 NANO

Automatic switching
among the 6 ranges

COM AUTO

Press any key on the keypad (except 1/0 (On/Off)) when the alphanumeric display indicates the desired range.

REMARK:

When the upper limit of the range is exceeded,

- 1) In locked range mode: the numerical display remains locked at 1999.
- 2) In automatic switching mode: the instrument switches to the next higher range.

After a power outage or being switched off and on, the C320 resets to the chosen range.

4.3.2 ADJUSTING THE CELL COEFFICIENT

We recommend adjusting the cell coefficient using a reference solution in the same range as the solution to be measured.

Immerse the conductivity cell and the temperature probe in the reference solution.

EXAMPLE: KCl 0.1 mol.L-1, range 0-20mS value 12.97mS at 25°C (see table: Conductivity of potassium chloride solutions, on page 40).

Press the:

key

The alphanumeric display alternates between:

COEF CEL

6 – 7 et 8

Press the:

or

key

To adjust the numerical display to the desired value

12.97

Press the:

key to validate.

The alphanumeric display indicates:

1.01 COEF CEL

English

The numerical display indicates the cell coefficient.
(The coefficient is expressed in cm: ratio of the surface area of the two plates of the cell to the distance between them).

The instrument then returns to the measurement mode.

4.3.3 VIEWING THE CELL COEFFICIENT

Press first the:

Then:


The numerical display indicates the cell coefficient:

1.01 **COEF CEL**

4.3.4 ADJUSTING THE TEMPERATURE COEFFICIENT

Because the solutions evolve differently as a function of temperature, it is necessary to set the temperature coefficient (expressed in % per degree) so as to refer the conductivity to 25°C whatever the actual measurement temperature.

Press the:

The alphanumeric display indicates:

COEF TEMP

Then:

6 - 7 et 8

The numerical display indicates the current value of this parameter (the factory setting is 2.2%/°C, the commonest value).

or

key.

Press the:

To set to the desired value.

Press on the:

key to validate.

The instrument then returns to the measurement mode.

NOTE: ADJUSTMENT RANGE: 0 to 8.

5. TABLE: CONDUCTIVITY OF POTASSIUM CHLORIDE SOLUTIONS

5.1 CONDUCTIVITY VS. TEMPERATURE

	Dilution to 1/10	Dilution to 1/50	Dilution to 1/100	Dilution to 1/1000
TEMPERATURE	KCl 0,1 mol.L ⁻¹	KCl 0,02 mol.L ⁻¹	KCl 0,01 mol.L ⁻¹	KCl 0,001 mol.L ⁻¹
°C	Conductivity mS/cm	Conductivity mS/cm	Conductivity μS/cm	Conductivity μS/cm
15	10.410	2.242	1 147	119.1
16	10.670	2.293	1 174	121.9
17	10.930	2.347	1 199	124.5
18	11.190	2.398	1 224	127.1
19	11.430	2.451	1 250	129.6
20	11.700	2.500	1 279	132.5
21	11.960	2.551	1 305	135.3
22	12.220	2.604	1 331	138.1
23	12.470	2.659	1 359	140.9
24	12.730	2.710	1 387	143.8
25	12.970	2.769	1 412	146.5

5.2 CONDUCTIVITY OF REFERENCE SOLUTIONS OF POTASSIUM CHLORIDE AT 20°C

Concentration (mol.L ⁻¹)	Dilution		Conductivity mS/cm	
1	Potassium chloride..... 74.56g Distilled water..... <i>q. s.p.</i> 1 000 mL		112.359	
	Volume of solution 1 mol.L ⁻¹ (mL)	Distilled water (mL)		
	0.400	40		<i>q. s.p.</i> 100
	0.200	20		100
	0.150	15		100
	0.100	10		100
	0.067	6.6		100
	0.050	5		100
	0.040	4		100
0.033	3.3	100	41.660	
0.025	2.5	100	22.070	
			16.390	
			11.700	
			7.812	
			5.988	
			4.807	
			4.032	
			3.048	
	Volume of solution 0,1 mol. L ⁻¹ (mL)	Distilled water (mL)	Conductivity µS/cm	
0.020	20	<i>q. s.p.</i> 100	2 510	
0.015	15	100	1 851	
0.010	10	100	1 279	
0.008	8	100	1 000	
0.0067	6.6	100	843.8	
0.005	5	100	653.5	
0.004	4	100	512.8	
0.0033	3.3	100	430.1	
0.0025	2.5	100	325.2	
0.002	2	100	298.5	
0.0015	1.5	100	196	
0.001	1	100	132.5	

Note: Take all essential precautions to obtain accurate measurements and homogeneous solutions.

5.3 ADJUSTING RECORDER OUTPUT

The conductivity recorder output is on the rear panel of the instrument:

- Black  and red  terminals for insulated 4mm banana plugs.

The recorder output is inoperative (0V) when the instrument is in the automatic range switching mode.

5.3.1 FACTORY OUTPUT SETTING

When first used or after a reset:

- 0% (0V) = 200
- 100% (5V) = 1800

5.3.2 ADJUSTING THE PARAMETERS OF THE ANALOG OUTPUT

Press the:

The alphanumeric display indicates in turn:

SORTIE

ANALOG

0/100 corresponds to 0V

0/100

The numerical display indicates the present setting.

Press the:

To adjust the numerical display to the desired 0/100 value (from 0 to 1999).

English

Then press on the:

key to validate

The alphanumeric display indicates:

100/100 corresponds to 5V

100/100

The numerical display indicates the present setting.

Press the:
validate.

or

key to

To adjust the numerical display to the desired 100/100 value (from 0/100 to 1999).

Then on the:

key to validate.

The C320 then returns to the measurement mode.


NOTE: Do not change the measurement range while using the recorder output.

5.4 RS 232/V24 SERIAL DIGITAL OUTPUT

This output makes data acquisition with a computer or a printer very simple.

CHARACTERISTICS OF THE SERIAL LINK

Link by 9-pin connector:

- Reception RXD ②
- Transmission TXD ③
- Request line ⑧
- Logical ground ⑤


Character format: 8 bits, no parity, 2 stop bits

Transmission speed: 2400 baud

Transfer of a string of 23 bytes, of which 10 bytes (9+CR) for the temperature and 13 bytes (12+CR) for the conductivity.

There are 2 transmission request modes:

1. Request by front: the command line must go from 5V to 0V.

2. Request by word: for each parameter, there is a code that can be used to select what is wanted.

Code required: T°C = E and Conductivity = C

An additional code is used to transfer the whole string:

Code required: G

5.5 RESET

Following a power outage or a magnetic or RF induction event, the memories may contain incorrect values.

The C 320 must then be reset to operate normally

5.5.1 RESET PROCEDURE:


Press the:

The alphanumeric display indicates:

VAL ****

Then the:

The alphanumeric display indicates:

VAL *123*
VAL *INIT

The C 320 returns to the measurement mode.

Following this reset, all of the adjustments must be redone.

5.5.2 VALUES AFTER RESET:

- Automatic range switching
- Cell coefficient: 1.00
- Coefficient %/°C: 2.2
- Manual temperature: 25.0°C
- Automatic temperature: no offset.
- Recorder output, 0/100 = 200
- Recorder output, 100/100 = 1800

6. IF SOMETHING GOES WRONG

PROBLEM	ADVICE
No display	<p>Check the 230 V supply.</p> <p>Check the connection of the power supply unit to the power connector</p>
No external temperature reading	<p>Check the connection of the temperature probe</p> <p>Check the value of the probe (107 Ω at 20 °C)</p>
Conductivity fails to stabilize in all solutions measured	<p>Check that the conductivity cell is correctly immersed.</p>
No variation of conductivity	<p>Check the connection of the cell.</p>
The instrument displays SOS-ALARM	<p>Check the condition.</p> <p>If necessary, replace it, or perform tests with pure resistances.</p> <p>For 1kΩ the value displayed 1000 μS.</p>
Displays 1999	<p>You are over the range limit.</p> <p>Switch to the next higher range.</p>

7. SAFETY RULES

- The mains supply must have the following characteristics: 230 V \pm 10% 50-60 Hz - 5W.
- The power supply unit serves as voltage disconnect device.
- The interior of the instrument must always be kept clean and dry.
- If the instrument is used in a way not in conformity with the specifications, the protection provided by the instrument may be impaired.
- Disconnect the instrument before opening the housing.
- The operation of the instrument may be perturbed by radiated electric fields or electrostatic discharges; operator intervention will be necessary to restore normal functioning.
In consequence, the instrument must not be used for permanent operation without human supervision.

8. MAINTENANCE

The instrument contains no parts that can be replaced other than by trained and accredited personnel. Any unauthorized repair or replacement of a part by an "equivalent" may gravely impair safety.

8.1 CLEANING

Use a cloth moistened with clean water or a neutral detergent to wipe the instrument, then wipe with a dry cloth.

Do not use the instrument again until it is completely dry.

8.2 METROLOGICAL CHECK

Like all measuring or testing devices, the instrument must be checked regularly.

This instrument should be checked at least once a year. For checking and calibration, contact one of our accredited metrology laboratories (information and contact details available on request), at our Chauvin Arnoux subsidiary or the branch in your country.

8.3 REPAIR

For all repairs before or after expiry of warranty, please return the device to your distributor.

9. TO ORDER

9.1 DELIVERY CONDITION

- 1 C 320 Conductimeter
- 1 user manual
- 125ml of 1Mol/l KCl solution
- 1 9V power supply unit

All packed in a cardboard box.

FRANCE

Chauvin Arnoux Group

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75876 PARIS Cedex 18

Tél : +33 1 44 85 44 85

Fax : +33 1 46 27 73 89

info@chauvin-arnoux.com

www.chauvin-arnoux.com

INTERNATIONAL

Chauvin Arnoux Group

Tél : +33 1 44 85 44 38

Fax : +33 1 46 27 95 69

Our international contacts

www.chauvin-arnoux.com/contacts